

Manual Crimping Tools



CRS's ergonomically curved handle design conforms to your natural hand grip, which greatly improves comfort over blocky metal grips. The body is constructed of tough, light-weight glassfilled plastic for strength.

The bottom-pull handle design is steadier to hold than top-push handle designs. The head does not wobble. It is easy to hold the tapered jaws of the crimper in a stable position, as the bottom handle is pulled up.

The adjustment knob is conveniently positioned at the head of the crimper, clearly displays + and - symbols along with directional arrows. This simplifies adjusting the intensity of the desired crimp. When the knob touches the crimper body, the crimp setting has been reached. No more "extra squeeze" necessary!

To ensure a smooth and consistent operation, the hand crimper uses the same precision tolerances as the CRS Electronic Crimper and Decapper.

Description	Item No.	Item No.
11mm Manual Crimper	311990	11mm Manual Decapper 311991
20mm Manual Crimper	320990	20mm Manual Decapper 320991
13mm Manual Crimper	313990	13mm Manual Decapper 313991
8mm Manual Crimper	308990	
		Base 5ASTA0
20mm Manual Crimper, for Flip Off Caps	320992	
13mm Manual Crimper, for Flip off Caps	313992	

Analytical Columns
Unit 7, Addington Business Centre
Vulcan Way, New Addington
Croydon. CR0 9UG



Tel : 01689 842736
Email : analytical@btconnect.com

www.analyticalcolumns.com